

# Certificate of Conformity

<b>VWR Catalogue Number</b>	See List
<b>Description</b>	Petri dishes, glass, VWR®, without vents, Non-Sterile
<b>Country of Origin</b>	Manufactured in India
<b>Date of Issue (yyyy-mm-dd)</b>	2022-05-30

Catalogue Number	Diameter [mm]	Height [mm]
391-0576	60	20
391-0577	80	20
391-0578	100	15
391-0579	100	20
391-0580	120	20
391-0581	150	25

## Quality System Compliance

Products are manufactured under the **ISO 9001 and ISO 14001** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

## QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

## Product Specifications

**Material** Type 1 Class A Borosilicate 3.3 glass.

**Sterilization** Non-Sterile.

**ATP Assay** Not Applicable.

<b>DNase &amp; RNase Free</b>	This product is free of any detectable DNase/RNase contamination.
<b>BSE/TSE</b>	Not Applicable.
<b>Cytotoxicity</b>	Not Applicable.
<b>Non-Pyrogenic Statement</b>	Not Applicable.
<b>Latex Statement</b>	Not Applicable.
<b>BPA Statement</b>	BPA is not intentionally used in the manufacture or the formulation of this product. However, this product has not been tested for this chemical material.
<b>DEHP Statement</b>	Not Applicable.
<b>RoHS</b>	Not Applicable.
<b>REACH Statement</b>	Not Applicable.
<b>Storage Conditions</b>	Room Temperature.
<b>Shelf Life</b>	2 years.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.  
This document has been produced electronically and is valid without a signature.